








Observation	Cq is Earlier than Expected	Variable Technical Replicates which have Cq Difference of >0.5 Cycles	Unexpected Data Values	No Data in Wells Selected for Analysis	Exponential Amplification of the No Template Control	Plateau is Lower than Expected	Comparisons between sample are Irreproducible
							
Potential Causes	<ul style="list-style-type: none"> The RNA has been contaminated by genome DNA There are multiple products There is a 'multicopy' gene Poor specificity of primer High primer-dimer production when using binding dye detection Transcript has high expression within the experimental samples 	<ul style="list-style-type: none"> Pipetting error Insufficient mixing of the solutions Low expression of target transcription resulting in stochastic amplification Reaction is poorly optimised There is a high Cq and low concentrations of the template 	<ul style="list-style-type: none"> Samples labelled incorrectly Plate has been inserted backwards Poor primer specificity Inhibitors in the sample 	<ul style="list-style-type: none"> Wells were not selected for analysis The wrong dye is selected for analysis First-strand synthesis failure PCR failure No expression of the target transcript 	<ul style="list-style-type: none"> There is contamination from laboratory exposure to the same target sequence There is contamination carried over from the reagent manufacture 	<ul style="list-style-type: none"> Reagents are limiting Reagents are degraded (dNTPs or master mix) Some probe dyes are less bright than others The reaction is inefficient Probe concentration is incorrect 	<ul style="list-style-type: none"> The amplification efficiency is below 88% in one or both samples Differences in efficiency are >5% Degradation of RNA Dilutions are inaccurate Low concentrations are being measured
Corrective Steps	<ul style="list-style-type: none"> The specificity of the primers require re-designing to be improved The assay performance needs to be tested carefully against quantified controls The primer concentration should be optimised The annealing temperature should be optimised DNase treatment should be used before reverse transcriptase 	<ul style="list-style-type: none"> Pipettes should be calibrated Use positive-displacement pipettes and filtered tips Mix all solutions thoroughly during preparation Hold pipette vertically when aspirating solutions - sterile technique does not ensure reproducibility when working with small volumes Optimise reaction conditions Add more sample to achieve lower Cq 	<ul style="list-style-type: none"> Rerun samples or plate with extra caution when loading Primers are re-designed to increase specificity Run a dilution of samples. Inhibitors will be diluted and may result in lower Cq for diluted material 	<ul style="list-style-type: none"> Settings for data collection and data viewing should be checked and reviewed Determine whether background fluorescence can be seen Repeat experiment with new reagents The assay performance should be tested against carefully quantified controls Check samples with an alternative target 	<ul style="list-style-type: none"> Work area should be cleaned with 10% bleach and nuclease-free water Prepare the reaction mix in a clean lab, separated from any template sources New reagent stocks should be ordered Ensure reagents are not contaminated by bacteria when using bacterial templates 	<ul style="list-style-type: none"> Calculations should be checked for master mix Repeat experiment using fresh stock solutions Compare end-point fluorescence to a different reaction using a probe labelled with the same dye 	<ul style="list-style-type: none"> Redesign primers for one or both genes Repeat experiment with fresh reagents and sample Check samples with a target that is more highly expressed